


Advanced X-ray Detector Technology




June 20, 2008



Confidential

1

The Electron Microscope



June 20, 2008




Confidential

2

The Electron Microscope

(The beam, the image formation and spectral x-ray detector)



The Beam

Focused beam of electrons.

Energy chosen between 1 and 30kV.

Electrons penetrate and lose energy in the specimen.

The Image Detector

Electrons ejected from sample are used for imaging in the SEI detector. (secondary electron detector)


The X-Ray Detector

X-rays emitted from sample are used for chemical analysis by the x-ray detector. The spectra is then analyzed for elemental identification

Information from Electron Microscope


Topography

The surface features of an object or "how it looks", its texture; direct relation between these features and materials properties




Morphology

The shape and size of the particles making up the object; direct relation between these structures and materials properties




Composition

The elements and compounds that the object is composed of and the relative amounts of them; direct relationship between composition and materials properties



Crystallographic Information

How the atoms are arranged in the object; direct relation between these arrangements and material properties



Evex X-ray Detectors Technologies



Evex - Si(Li)
Silicon Lithium
Crystal
129 FWHM
Lower Limit Be Beryllium
Window Size 10, 30, 50, 100mm



Evex – Ge(Li)
Germanium Lithium
Crystal
125 FWHM
Lower Limit Be (Beryllium)
Window Size 10, 30, 50, 100 mm



Evex-QDD Violin
Silicon
Chip Based
130 FWHM
Lower Limit Be – Beryllium
Window Size 5, 20, 50, 100 mm

X-ray Detector Windows - Thresholds

	SSUTW Super (UTW)	SUTW	UTW Ultra Thin	TW or (Be) Beryllium
Lower Detection Element	Be (4)	B (5)	C (6)	Na (11)
Energy	.109	.185	.282	1.041




SSUTW – Super Super Ultra Thin Window

SUTW – Super Ultra Thin Window

UTW – Ultra Thin Window

TW or (Be) – Beryllium Window



X-ray Windows

X-Ray Transmission

Light Rejection

Vacuum Tightness



Evex SSUTW windows are the highest performing x-ray windows available for low energy x-ray analysis.

Evex UTW windows are ideal for applications that require high x-ray transmission of light element energies, high mechanical strength, light rejection, vacuum tightness, and reliability.

June 20, 2008



Confidential

7

Si(Li)

vs

QDD



Evex - Si(Li)


- Detection
 - Sodium and Above
 - Beryllium and Above
- Windows
 - Beryllium (Be)
 - Light Element (UTW)
 - Super Light Element (SUTW)
- Count Rate <10,000 cps
- Peak Stability at Any Count Rate
- 129 FWHM
- Liquid Nitrogen Cooling
- 2-4 of cooling before use



Evex-QDD Violin

- Detection
 - Sodium and Above
 - Boron and Above
- Windows
 - Beryllium (W) – Lowest Sodium - Na
 - Light Element (UTW) – Carbon - C
 - Super Light Element (SUTW) – Boron - B
 - Super Super Light Element (SSUTW) – Beryllium - Be
- Count Rate > 250,000cps
- Peak Stability at any Count Rate
- 130 FWHM
- No Liquid Nitrogen required for cooling
- Instant On

The Ultimate X-ray NanoAnalysis Machine For Electron Microscopy



Advanced Detector Technology





- High Count Rate
- Light element detection
- No liquid nitrogen
- No vibrations
- No maintenance
- Great for high magnification
- Compact design

Benefits of Evex-QDD

- **High Count Rate**
- **Light element detection**
- **No liquid nitrogen**
- **No vibrations**
- **Low maintenance**
- **Great for high magnification**
- **Compact design**
- **Trace Analysis**



Inside the Evex QDD – Liquid Nitrogen Free, High Speed, X-ray Detector



June 20, 2008



Confidential

11

Evex Mini-SEM with X-ray NanoAnalysis




Website
<http://www.mini-sem.com>

Electron Microscope Specification

Magnification	30,000 x (Digital Zoom x 4 ~ 120,000 x)
Acceleration Voltage	1 kv ~ 30 kv (variable)
Detectors	SEI, BEI, EDS
X-ray Detector	Evex QDD – Liquid Nitrogen Free, High Speed, X-ray Detector
Image Size	4096 x 4096
Image Format	JPG, TIFF, BMP, PNG
Data Display	Micron Bar, Magnification. Date, Kv
Voltage	110 volts / 220 volts

Evex NanoAnalysis System



June 20, 2008




Confidential

14

Workstation

Workstation

3.0 GHZ Processor
1.0 GB Memory
Microsoft Windows XP OS
160 GB Hard Drive
22" LCD Display
Network Card
Microsoft Office - Optional



June 20, 2008



Confidential

15

NanoAnalysis Features

Image Acquisition Module

- Image Processing
- Image Analysis
- Image Operations
- Image Filters
- Thermal Imaging
- Image Measurement / Annotation
- Easy Report Writing
- AutoSave




Spectral Acquisition Module

- Spectral Processing
- Spectral Analysis
- Qualitative Analysis
- Quantitative Analysis
- Monte Carlo Simulation
- Detector Geometry Wizard
- Auto ID with Peak Labeling
- Auto Region of Interest (ROI)
- Automatic Calibration & Logging
- Easy Report Writing
- AutoSave

Elemental Mapping

- Standard Elemental Mapping
- Ultra Fast Mapping
- Line Scan
- Point Mode & MultiPoint
- AutoAnalysis (Spectra at Every Pixel)
- ColorSEM - Quick Rendering
- Focus Mapping
- Thermal Mapping
- Easy Report Writing
- AutoSave



Advanced Features

- Chemical Typing
- Trace Sensitive Analysis
- Critical Measurement
- 3D Mapping
- 3D Measurement
- AutoAnalysis (Spectra at Every Pixel)
- Remote Conferencing
- Stereo Imaging
- Synchro-Analysis