

# dPette+

## Multifunctional 8-channel Electronic Pipette



DLAB multifunctional eight-channel electronic pipette is an electric-driven pipetting device with digital controls. Our latest electronic technology assures high efficiency, accuracy and precise results for modern regulated laboratory reducing human error regardless of user skill level, and is a better choice for scientific researchers.



### Function Wheel

Switch functions such as pipetting, continuous dispensing, and pipetting speed adjustment

### Parameter Knob

Parameter adjustment  
Parameter confirm  
Aspirate/ Dispense/Mix

### Motor Driven Digital Control

### 360° Rotation

### Ingenious Curved Ejector

Tips are pushed out instantaneously from outside to inside to minimize the thrust.

### ● High accuracy and precision, time-saving and efficient

High quality stepper motor and piston control system provide outstanding accuracy and repeatability results, which reduce manual operation errors.

### ● High level ergonomic and comfortable handle design

Lightweight and compact, to protect the user from repetitive strain injury (RSI).

### ● Simple to use and versatile

Double knobs for simple and versatile control. Dispensing function reduces the pipetting steps and speeds up work.

### ● Dual charging mode

Using the USB charger or charging stand to ensure uninterrupted use.

### Three steps to operate dPette+ easy



Long press the **parameter knob** for 2s to start




Quickly turn the **function wheel** to activate, rotate to switch pipetting, continuous dispensing, and other function settings



After quickly turning the **parameter knob** to unlock, turn to adjust the volume, press to pipette, long press to mix


## Pipetting function

Parameter knob press → Pipetting,  
Long press → Mixing




## Dispensing function

Dispensing times = Pipette max volume /  
Single sample volume



## Other function

Pipetting speed adjustment  
Key tone adjustment



## Complimentary pipette holder

Supports dual charging mode



### USB charging



Support automatic calibration  
Connect to the computer,  
with free software provided

### Contact charging



## Specifications

Volume Range	Increment	Test Volume	Accuracy error		Precision error	
µL	µL	µL	µL	%	SD* µL	CV* %
0.5-10 µL	0.01	10	±0.20	±2.00	0.1	1.0
		1	±0.08	±8.0	0.05	5.00
10-100 µL	0.1	100	±0.80	±0.80	0.3	0.30
		10	±0.3	±3.00	0.2	2.00
30-300 µL	1	300	±1.80	±0.60	0.9	0.3
		30	±0.90	±3.00	0.3	1.0

\* SD=Standard Deviation \* CV=Coefficient of Variation

\*DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655.

