

NMR Tubes, **5 mm**, Disposable Grade

Article Number	Wilmad Equivalent	Approx. MHz	Length, Inches	Wall Thickness
897193-0000	512-7	60	7	0.015
897193-0008	512-8	60	8	0.015
O.D. INCHES: 0.193-0.196			I.D. INCHES: 0.163-0.168	
897194-0000	512-7	60	7	0.02
O.D. INCHES: 0.190-0.197			I.D. INCHES: 0.150-0.157	

NMR Tubes, **3 mm**,
ASTM Type 1, Class A, Borosilicate Glass (KG-33)

Article Number	Wilmad Equivalent	Approx. MHz	Length, Inches	Wall Thickness	Wall Variation T.I.R.†	Camber
897840-0000	335-PP-7	>400	7	0.01118	0.0005	0.00025
897840-0008	335-PP-8		8	0.01118	0.0005	0.00035
897840-0009	335-PP-9		9	0.01118	0.0005	0.00045
897835-0000	328-PP-7	360	7	0.01118	0.0010	0.0005
897835-0008	328-PP-8		8	0.01118	0.0010	0.0006
897835-0009	328-PP-9		9	0.01118	0.0010	0.0007
897830-0000	327-PP-7	270	7	0.01118	0.0015	0.0010
897830-0008	327-PP-8		8	0.01118	0.0015	0.0011
897830-0009	327-PP-9		9	0.01118	0.0015	0.0012
897825-0000	—	200	7	0.01118	0.0020	0.0005
897825-0008	—		8	0.01118	0.0020	0.0006
897825-0009	—		9	0.01118	0.0020	0.0007
897820-0000	307-PP-7	150	7	0.01118	0.0020	0.0010
897820-0008	307-PP-8		8	0.01118	0.0020	0.0011
897820-0009	307-PP-9		9	0.01118	0.0020	0.0012
897805-0000	—	100	7	0.01118	0.0025	0.0020
897805-0008	—		8	0.01118	0.0025	0.0020
897805-0009	—		9	0.01118	0.0025	0.0020
897800-0000	305-PS-7	80	7	0.01118	0.0030	0.002
897800-0008	305-PS-8		8	0.01118	0.0030	0.002
897800-0009	305-PS-9		9	0.01118	0.0030	0.002

O.D. INCHES - 0.1186" (+0.0000/-0.0005) I.D. INCHES - 0.0945" (+0.0005/-0.0000) † T.I.R. = TOTAL INDICATOR READING IN INCHES.

NMR Tubes, **10 mm**,
ASTM Type 1, Class A, Borosilicate Glass (KG-33)

Article Number	Wilmad Equivalent	Approx. MHz	Length, Inches	Wall Thickness	Wall Variation T.I.R.†	Camber
897335-0000	513-7PP-7	>360	7	0.0182	0.0015	0.0005
897330-0000	513-5PP-7	150	7	0.0182	0.0020	0.0010
897325-0000	513-3PP-7	80	7	0.0182	0.0030	0.0015
897320-0000	513-1PP-7	60	7	0.0182	0.0090	0.0020

O.D. INCHES - 0.3937 (+0.0000/-0.0005) I.D. INCHES - 0.3569 (+0.0005/-0.0000) † T.I.R. = TOTAL INDICATOR READING IN INCHES.