

## NMR Tube Caps, PTFE, 3 and 5 mm

NMR Tube caps for use with any Kontes 5 mm standard or Medallion series NMR tube. These PTFE caps are precision machined and offer the user a means of universal solvent compatibility. The caps are sold in packages of 25.

Article Number	Fits Tube Size (mm)	Color	Pkg. Qty.
897096-0003	3	White	25
897096-0005	5	White	25

## NMR Tube Pressure Caps, PE, 3, 5 and 10 mm

A tight-fitting, polyethylene cap which stays in place even with a build-up of internal vaporization pressures.

Article Number	Wilmad Equivalent	Fits Tube Size (mm)	Color	Pkg. Qty.
897093-0001	—	3	Red	100
897095-0001	521	5	Red	100
897095-0021	521	5	Black	100
897095-0061	521	5	White	100
897095-0071	521	5	Green	100
897095-0081	521	5	Yellow	100
897095-0011	—	5	Assorted	100
897095-1001	521-C	10	Red	100
897095-1021	—	10	Black	100
897095-1071	521-C	10	Green	100
897095-1081	521-C	10	Yellow	100
897095-1011	—	10	Assorted	100

## NMR Tube Septum, Rubber

Useful as a closure or for sample outgassing. Soft, natural rubber allows easy penetration with a syringe needle.

Article Number	Wilmad Equivalent	Fits Tube Size (mm)	Pkg. Qty.
897097-0001	521-S	5	100

## NMR Tube Sample Filter

Pressure-type filter for removal of particulate matter in samples prior to analysis. Upper part is fitted with a #00 stopper. Sample is pipetted into the filter. A rubber bulb, with a valve, provides pressure to force the sample through the fine porosity filter. 897050-0021 has a bottom with a plain drip tip. 897050-0022 has an inner luer joint as a tip.

Article Number	Wilmad Equivalent	Style	Tip Configuration	Pkg. Qty.
897050-0021	807	21	Plain Drip Tip	1
897050-0022	808	22	Inner Luer Drip Tip	1

## NMR Tube Funnel

Polypropylene construction with a lower fitting sized to accommodate 5 mm NMR tubes. Disposable feature eliminates the possibility of cross contamination. Filter is a polyethylene disc or glass wool. Complete set consists of a reservoir, a column and a lower fitting.

Article Number	Description	Pkg. Qty.
420160-0000	Funnel Set (Unassembled)	100

### Replacement Parts

420164-1000	Reservoir Only	50
420168-1000	Lower Fitting Only	50

### Accessories

420162-0000	Polyethylene Disc, 30-50 Microns	100
420162-0020	Polyethylene Disc, 20 Microns	100