

# LARGE VOLUME NMR SAMPLE TUBES

## 8mm NMR SAMPLE TUBES (0.3145 - 0.3150" OD) ROUND BOTTOM

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-24	25+
NE-L8-7 NE-L8-8	100-300	0.50mm	0.003"	0.002"	7"(178mm) 8"(203mm)	18,05 19,35	17,05 18,35
NE-H8-7 NE-H8-8	400 -500	0.50mm	0.0015"	0.001"	7"(178mm) 8"(203mm)	29,25 31,15	27,90 29,65

## 10mm NMR SAMPLE TUBES (0.3932 - 0.3937" OD) ROUND BOTTOM

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-24	25+
NE-L10-7 NE-L10-8	100 - 300	0.50mm	0.003"	0.002"	7"(178mm) 8"(203mm)	14,15 15,70	13,55 15,15
NE-H10-7 NE-H10-8	400 -500	0.46mm	0.0015"	0.001"	7"(178mm) 8"(203mm)	23,50 25,65	23,00 24,55
NE-U10-7 NE-U10-8	600+	0.46mm	0.001"	0.0005"	7"(178mm) 8"(203mm)	28,05 29,45	25,50 26,80
NE-H10-MW-7	400	1.00mm	0.0015"	0.001"	7"(178mm)	27,45	24,95
NE-H10-HW-7	400	1.7mm	0.003"	0.001"	7"(178mm)	32,90	29,65

# HIGH TEMPERATURE / UV STUDIES

## 10mm QUARTZ NMR SAMPLE TUBES (0.3932 - 0.3937" OD) ROUND BOTTOM

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-24	25+
NE-MQ10-7	200	0.50mm	0.003"	0.002"	7"(178mm)	79,65	75,65
NE-HQ10-7	300 -500	0.50mm	0.0015"	0.001"	7"(178mm)	81,30	77,25

## 15mm NMR SAMPLE TUBES (0.5900 - 0.5910" OD) ROUND BOTTOM

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-9	10+
NE-L15-7	200-400	0.70mm	0.002"	0.001"	7"(178mm)	32,45	28,75
NE-L15-M-7		1.30mm	0.003"			Request	Request

## 16mm NMR SAMPLE TUBES (0.6294 - 0.6304" OD) ROUND BOTTOM

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-9	10+
NE-L16-7	200-400	0.72mm	0.002"	0.001"	7"(178mm)	Request	Request
NE-L16-8					8"(203mm)	Request	Request

Note: 20-30mm OD sample tubes are suitable for imaging applications

Lengths and Bottom Style are optional.

**20mm NMR SAMPLE TUBES (0.7870 - 0.7880" OD) FLAT BOTTOM**

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-9	10+
NE-L20-7	200 - 400	0.80mm	0.002"	0.001"	7"(178mm)	48,81	46,25
NE-L20-8					8"(203mm)	54,35	51,65

**25mm NMR SAMPLE TUBES (0.9838 - 0.9848" OD) FLAT BOTTOM**

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-9	10+
NE-L25-7	200-400	0.90mm	0.002"	0.001"	7"(178mm)	58,40	55,80
NE-L25-8					8"(203mm)	62,05	58,15

**27mm NMR SAMPLE TUBES (1.0625 - 1.0635" OD) FLAT BOTTOM**

CATALOG NUMBER	Usage MHz	Wall	Maximum Wall Variation	Camber TIR	Length	Price € / each	
						1-9	10+
NE-L27-7	200-400	0.90mm	0.002"	0.001"	7"(178mm)	62,05	58,15
NE-L27-8					8"(203mm)	67,45	64,10

**30mm NMR SAMPLE TUBES (1.18 ± 0.02" OD) ROUND OR FLAT BOTTOM**

CATALOG NUMBER	Usage	OD	Wall	Length	Price € / each
NE-SL30	Imaging	30±0.5mm	1.4mm	optional	Request