

NMR SAMPLING ACCESSORIES

Catalog Number	Tube OD	Description	Price €/each
NE-230-3	3mm	straight body, 10mm OD	57,00
NE-230-3J	3mm	with 24/40 inner joint	70,35
NE-230-5	5mm	straight body, 10mm OD	57,00
NE-230-5J	5mm	with 24/40 inner joint	70,35
NE-230-10	10mm	straight body, 14mm OD	68,10
NE-230-10J	10mm	with 24/40 inner joint	69,30

NMR SAMPLE TUBE WASHERS

VACUUM FLASK

Catalog Number	Description	Price €/each
NE-240	Vacuum Flask, 500ml	60,20

MULTIPLE SAMPLE TUBE WASHER, 3 and 5mm versions

Catalog Number	Description	Price €/each
NE-232-3	3mm - 5 Tube Washer, complete	692,90
NE-232-5	5mm - 5 Tube Washer, complete	609,40

NMR SAMPLE TUBE WASHERS

Catalog Number	Description	Price €/each
NE-231-RC3	3mm NMR Tube Washer, Teflon	142,85
NE-231-RC5	5mm NMR Tube Washer, Teflon	145,25
NE-231-RC10	10mm NMR Tube Washer, Teflon	162,70
NE-240	Vacuum Flask, 500ml	60,20

Standard Components

Catalog Number	Description	Price €/each
NE-231-B	Base only, PTFE Reservoir only, PTFE	73,85
NE-231-R		44,80

3mm Components

Catalog Number	Description	Price €/each
NE-231-C3	Cap Gasket, 3mm, PE	2,00
NE-231-G	Glass Body only	10,05
NE-231-S3	Solvent Capillary, PEEK, 1/16"	7,35
NE-231-N3	Nut, PPS, 1/16"	7,70
NE-231-F3	Ferrule, ETFE, 1/16"	4,15

5mm Components

Catalog Number	Description	Price €/each
NE-231-C5	Cap Gasket, 5mm, PE	2,00
NE-231-G	Glass Body only	10,05
NE-231-S5	Solvent Capillary, PEEK, 1/8"	11,30
NE-231-N5	Nut, PPS, 1/8"	7,70
NE-231-F5	Ferrule, ETFE, 1/8"	4,15

10mm Components

Catalog Number	Description	Price €/each
NE-231-C10	Cap Gasket, 10mm, PE	2,85
NE-231-G10	Glass Body only	33,50
NE-231-S5	Solvent Capillary, PEEK, 1/8"	11,30
NE-231-N5	Nut, PPS, 1/8"	7,70
NE-231-F5	Ferrule, ETFE, 1/8"	4,15

SAMPLE TUBE CAPS

Catalog Number	Tube Size	Description	Price	
			Pk/100	Pk/1000
NE-310-3	3mm	LD Polyethylene, Red only	7,60	58,60
NE-310-5	5mm	EVA; Blue, Green, Red, Yellow, White, Black; specify color	10,90	72,45
NE-310-5S	5mm	Serum Stopper, sleeve type	56,60	-
NE-310-8	8mm	LD Polyethylene, Neutral only	15,80	92,75
NE-310-10	10mm	LD Polyethylene, Red only	15,80	92,75
NE-310-15	15mm	LD Polyethylene, Red only	39,05	339,15
NE-310-16	16mm	LD Polyethylene, Red only	39,05	339,15
NE-310-20	20mm	LD Polyethylene, Neutral only	39,05	339,15
NE-310-25	25mm	LD Polyethylene, Red only	39,05	339,15
NE-310-27	27mm	LD Polyethylene, Red only	39,05	339,15

5mm TEFLON SAMPLE TUBE CAPS

Catalog Number	Description	Price €
NE-312-5-10	Cap, PTFE, 0.8mm wall x 15mm long with circular external groove	54,00/Pk/10
NE-312-5-25		124,15/Pk/25
NE-312-5-50		212,95/Pk/50

PIPETS / RUBBER BULBS

Catalog Number	Description	Price €
NE-201	Long Tip Pipet for 5mm x 7" tubes; 7.5" tip, 10.5" OAL, Bx/100	98,20
NE-201-A	Extra Long Tip Pipet for 5mm x 8-9" tubes; 9.25" tip, 12.5" OAL, Bx/100	102,15
NE-201-B	Long Tip Pipet for 3mm x 7" tubes; gaged to fit to the bottom of the tube; 7.5" tip 10.5" OAL, Bx/100	102,15
NE-301-50	Rubber Bulbs, Pk/50	37,55

SAMPLE TUBE HOLDER

Catalog Number	Size	Description	Price €/each	
			1-4	5+
NE-330-5	8" x 4" x 9"	Sample Tube Holder, 5mm x 72 positions	53,70	51,40

CLEANING BRUSH

Catalog Number	Description	Price €/each	
		1-4	5+
NE-343-5	Brush, for 5mm OD thin-wall tubes, 50mm of bristles, 240mm OAL	11,00	10,40

MICRO SPATULA

Catalog Number	Description	Price €/each
NE-342	Micro Spatula - Stainless Steel	46,00

SAMPLE TUBE LABELS

Catalog Number	Description	Price €/roll
NE-7407	NMR Tube Labels, Clear, for 5mm OD Tubes	48,00

MARKING PENS

Catalog Number	Description	Price €/each
NE-7301	Blue ink, fine point, fast drying	4,15
NE-7302	Black ink, fine point, fast drying	4,15
NE-7303	Red ink, fine point, fast drying	4,15
NE-7304	Green ink, fine point, fast drying	4,15

QUARTZ PROBE INSERTS

Catalog Number	Probe Size	Description	Price €/each
NE-295-5	5mm	Precision, thin wall quartz sleeve, 5.8mm x 6.8mm x 64mm. Replacement for Bruker spectrometers. Conc.= 0.001" TIR.	139,25
NE-295-5A	5mm	Precision, thin wall quartz sleeve, 5.6mm x 6.1mm x 64mm. For late model Bruker probes. Conc.= 0.001" TIR.	98,55
NE-295-5B	5mm	Precision, thin wall quartz sleeve, 5.6mm x 6.1mm x 54.5mm. For Bruker dual flow probes. Conc.= 0.001" TIR.	98,55
NE-295-V-5	5mm	Precision, thin wall quartz sleeve, 5.5mm x 6.0mm x 58.6mm. For Varian probes. Conc.= 0.001" TIR.	139,25
NE-295-12	12mm	Precision, thin wall quartz sleeve, 11mm x 12mm x 75mm. For micro imaging systems.	61,80
NE-295-22.8	22.8mm	Precision, thin wall quartz sleeve, 21mm x 22.8mm x 75mm. For micro imaging systems.	71,60
NE-295-30	30mm	Precision, thin wall quartz sleeve, 28mm x 30mm x 75mm. Replacement for Bruker micro imaging spectrometers.	96,20