

NMR MICRO TUBES

Micro NMR Sample Tubes

CATALOG NO.	APPLICATIONS	PRICE (€)/EACH		
		1-4	5-9	10+
NE-H5/4	Varian 5mm probes w/wo auto-changer. Bruker and JEOL 5mm probes.	37,00	34,00	32,35
NE-H5/4-Br	Bruker 5mm probes w/wo auto-changer. JEOL and Varian 5mm probes	37,00	34,00	32,35
NE-H5/3	Varian 3-5mm probes w/wo auto-changer. Bruker and JEOL 5mm probes.	37,00	34,00	32,35
NE-H5/3-Br	Bruker 5mm probes w/wo auto-changer. JEOL and Varian 5mm probes.	37,00	34,00	32,35
NE-H5/2.5	Bruker 2.5-5mm probes w/wo auto-changer. Varian 3-5mm probes w/wo auto-changer.	37,00	34,00	32,35
NE-H5/2.5-V	Bruker and JEOL 5mm probes. Bruker 2.5-5mm probes w/wo auto-changer.	37,00	34,00	32,35
NE-H5/2.5-V-Br	JEOL 5mm and Varian 5mm probes w/wo auto-changer.	37,00	34,00	32,35
NE-H5/2-V-Br	Bruker 2.5-5mm probes w/wo auto-changer. JEOL 5mm and Varian 3-5mm probes.	37,00	34,00	32,35

All volumes are calculated at 30 mm sample height and are for reference, only.
Proper filling factors may vary with the probe.

Capillaries

Capillary Size Mm	Volume At 30mm Height	Cat.No. 75mm Long	PRICE (€)/PK/10			Cat.No 100mm Long	PRICE (€)/PK/10		
			1-4	5-9	10+		1-4	5-9	10+
1 x 0.8	15µl	NE-262-1	47,25	43,95	41,80	NE-263-1	58,95	54,45	51,55
1 x 0.8-PS	15µl	NE-262-1-PS	64,70	60,20	57,20	NE-263-1-PS	75,15	71,50	67,85
1.7 x 1.4	50µl	NE-262-1.7	47,25	43,95	41,80	NE-263-1.7	58,95	54,40	50,95
2 x 1.2	35µl	NE-262-2	47,25	43,95	41,80	NE-263-2	58,95	54,40	50,95
2.5 x 2	97µl	NE-262-2.5-A	47,25	43,95	41,80	NE-263-2.5A	58,95	54,40	50,95
2.5 x 2.2	115µl	NE-262-2.5	47,25	43,95	41,80	NE-263-2.5	58,95	54,40	50,95

ADAPTERS

CATALOG NUMBER	FOR USE WITH THESE CAPILLARIES	Price (€)/each		
		1-4	5-9	10+
NE-325-5/1	NE-262-1-PS, NE-263-1-PS	36,20	35,60	35,00
NE-325-5/1.7	NE-262-1.7, NE-263-1.7	36,20	35,60	35,00
NE-325-5/2	NE-262-2, NE-263-2	36,20	35,60	35,00
NE-325-5/2.5	NE-262-2.5, -2.5A; NE-263-2.5, -2.5A	36,20	35,60	35,00

EXTRACTION ROD

CATALOG NUMBER	DESCRIPTION	Price (€)/each		
		1-4	5-9	10+
NE-341-5	Extraction Rod, Kel-F, 1/8 X 8", 4-40 thread	23,25	22,15	21,05

MICRO TUBE / CAPILLARY KITS

Each kit contains one micro tube/sample tube, one adapter, one extraction rod and 5 capillaries

CATALOG NUMBER (Kit)	Micro Tube	Capillary	Capillary OD	Capillary Volume	Annular Volume	Price (€)/each
NE-380-A	NE-H5/3	NE-262-2	2mm	35µl	36µl	103,25
NE-380-B	NE-H5/2.5	NE-262-1.7	1.7mm	50µl	47µl	103,25
NE-380-C	NE-UP5-7	NE-262-2.5	2.5mm	115µl	268µl	85,15
NE-380-D	NE-H5/3-Br	NE-262-2	2mm	35µl	36µl	103,25
NE-380-E	NE-H5/2.5-V-Br	NE-262-1.7	1.7mm	50µl	47µl	103,25

All volumes are calculated at **30 mm sample height** and are for reference, only.
Proper filling factors may vary with the probe.

AUTO PREP / SAMPLING SYSTEMS - SAMPLE TUBES / MICRO

Catalog Number	Sample Volume (µl)	Body OD	Stem OD	Price (€)/each		
				1-24	25-99	100+
NE-HL5-4-AS	415	5mm	-	Inquire	Inquire	Inquire
NE-H5/2-AS	35	5mm	2.0mm	19,55	18,50	17,60
NE-H5/2.5-AS	115	5mm	2.5mm	19,55	18,50	17,60
NE-H5/3-AS	135	5mm	3.0mm	19,55	18,50	17,60

OFF-COIL LOCK

Catalog Number	Lock Volume	Price (€)/each		
		1-9	10+	
NE-5-OCL	~110 µl	30,30	28,65	
NE-5-CIC	~60 µl	29,25	28,65	
		1-4	5-9	10+
NE-H5/3	Outer Tube	37,00	34,00	32,35

1.7mm MICRO NMR CAPILLARY

Catalog Number	Description	Price (€)/Vial/25	
		1-4	5+
NE-1.7-125	Micro NMR Capillary, 1.7mm OD	44,70	40,25

3mm NMR SAMPLE TUBES (0.1175" - 0.1180")

Catalog Number	Description	Length	Price (€)/each	
			1-9	10+
NE-H3-7	400+ MHz; 3mm OD x 2.36mm ID x 0.32mm Wall	7" (178mm)	10,30	9,85
NE-H3-8	Filling Factor: 118 µl per 27mm height Concentricity 0.001" TIR, Camber 0.001" TIR	8" (203mm)	11,65	11,30